**SAMSUNG**

**SAMSUNG AMERICA**

Made in Taiwan

105 Challenger Road, 6th Floor  
Ridgefield Park, NJ 07660

[www. SAMUNGHealthyLiving.com](http://www.SAMUNGHealthyLiving.com)

AUTOMATIC

## Automatic Inflate Blood Pressure Monitor

- Large Three Row Digital Display
- Memory Stores 60 Measurements
  - Time and Date Indications
  - Pulse Display
  - AC Capability

## **Important Warnings/Medical Disclaimer**

Before using your blood pressure monitor, you should read and understand all instructions and follow all warnings.

Information in this manual is provided for informational purposes only. This manual and product are not meant to be a substitute for the advice provided by your own physician or other medical professional. You should not use the information contained herein or this product for diagnosing or treating a health problem, disease, or prescribing any medication. If you have or suspect that you have a medical problem, promptly contact your healthcare provider.

## In This Manual

Healthy Living®	3-4
Commonly Asked Questions About Blood Pressure	5-8
Intended Use	8
Measurement Method	9
Main Unit and Accessories	10
Special Features	11
Accuracy	12
Special Conditions	12
Digital Display Symbols	13
Battery Installation	13
AC Capability	14
Selecting Time and Date	14
Important Reminders	15-16
Applying the Cuff	17
Taking a Blood Pressure Measurement	18-19
Sixty Measurement Memory	20
Taking Care of Your Blood Pressure Monitor	21
Troubleshooting	22-23
Five Year Limited Warranty	24-25
Specifications	26
Blood Pressure Diary	27-30

## Healthy Living®

This blood pressure monitor is the first step to living a healthier life. But monitoring your blood pressure at home is just the beginning, it's also important to:

**EXERCISE REGULARLY** - Start out by taking a quick walk each day. Before you know it, you'll feel energized. That's because along with being an easy way to shed a few pounds, walking also helps bones and heart to become stronger.

**EAT RIGHT** - Remember, you are what you eat! So when you reach for a snack, try grabbing a piece of fruit or a vegetable. When cooking meals at home try to use fats and oils sparingly. Sometimes a healthy diet is the best medicine of all! It's also a good idea to monitor your sodium intake since about one-third of individuals who have high blood pressure are also sodium sensitive. Consult your doctor regarding over-the-counter supplements and herbs such as onion and garlic. Both of these plants have been historically linked to producing mild anti-hypertension effects.

## Healthy Living<sup>®</sup>

**BE HAPPY!** - A positive attitude is a big part of a healthy body. Try not to get too stressed out over things. Daily relaxation techniques such as gentle yoga or even quiet time by yourself can do a lot of good for the mind and body. Stay focused on the good things in life. Begin by congratulating yourself for that short walk you just took or the smart food choice you made. It's your first step to a healthier life!

**KEEP AN EYE ON YOUR WEIGHT!** - Surprisingly, losing even a small amount of weight can often times help control high blood pressure. Talk with your doctor to find a "healthy" weight for you and strive to come within a range of it.

**Remember** to consult your physician prior to beginning exercise, vitamin/mineral, dietary changes or programs and blood pressure monitoring.

## Commonly Asked Questions About Blood Pressure


### ***Why Should I Monitor My Blood Pressure At Home?***

By monitoring at home, you can manage your health between doctor visits. Home monitoring also provides you with a record of blood pressure measurements over time to give to your doctor - a valuable tool in monitoring hypertension.


Home monitoring may also help alleviate *White Coat Hypertension*. This term refers to the medical phenomena which causes many individuals' blood pressure to be higher when taken in a medical office setting. That's because many times the stress of being in the doctor's office may cause blood pressure to rise as much as 20-50 points higher than normal. When monitoring at home, you may get a truer indication of your blood pressure in a normal, relaxed state.

### ***What Is Blood Pressure?***

The commonly used term "blood pressure" refers to the force or *pressure* that is exerted on the body's arteries as blood flows through them. Each time the heart beats or *contracts*, it produces pressure in the arteries. This blood pressure is what moves blood through the body, supplying oxygen and nutrients to every organ.


## Commonly Asked Questions About Blood Pressure

When the heart contracts it sends blood out into the body. This pressure is referred to as **systolic**.

The pressure that is created as the heart pulls blood back through is called **diastolic**.


Blood pressure is recorded as systolic over diastolic, as shown here.

### ***What Is Considered Normal Blood Pressure?***

The Seventh Report of the Joint National Committee on Prevention, Detection, Evaluation, and Treatment of High Blood Pressure from the National Heart, Lung, and Blood Institute has set forth the following guidelines regarding high blood pressure. Guidelines provided are for information purposes only. Please consult with your physician for proper diagnosis.

## Commonly Asked Questions About Blood Pressure

### ***Why Does My Blood Pressure Fluctuate?***

Because blood pressure changes with every beat of the heart it is in constant fluctuation throughout any given day. In addition to these natural changes there are other factors that may affect blood pressure such as:

- Temperature (too hot or cold)
- Humidity
- Restlessness
- Lack of sleep
- Physical exercise
- Constipation
- Stress
- Fear
- Anger
- Anxiety
- Food and beverage consumption
- Individual physical conditions

In addition to these factors you should talk with your physician to see what may be affecting you.

### ***How Do I Know I Have High Blood Pressure?***

It's nearly impossible to know if your blood pressure is elevated without having it checked. People have the common misconception that one can "feel" their blood pressure rise. It is because of

## Commonly Asked Questions About Blood Pressure

this mystery that hypertension has been named “the silent killer.” Don’t let it sneak up on you! Start monitoring your blood pressure every day with your new Healthy Living<sup>®</sup> monitor.

### ***Why Is It Important To Control My Blood Pressure?***

Left uncontrolled, hypertension can place an incredible amount of stress on your heart and arteries, both of which become forced to work harder to keep blood flowing in a healthy manner throughout the body. Over time, this stress may result in health problems including over-enlargement of the heart, loss of elasticity of the arteries, heart disease, stroke and even death. Making some simple lifestyle changes now may help you to stay healthy in the future.

## Intended Use

Healthy Living<sup>®</sup> blood pressure monitors are convenient for home use. They are noninvasive, meaning no part of the monitor enters your body, and provide systolic, diastolic and pulse measurements. Healthy Living<sup>®</sup> monitors are recommended for use by people over the age of 18 and are not suitable for clinical use.

## Measurement Method


Healthy Living<sup>®</sup> monitors measure blood pressure and pulse by using what’s called an *oscillometric* method, meaning they measure the fluctuations in pressure. This is done by first wrapping a fitted cuff around the upper arm. The monitor then automatically fills the cuff with air creating pressure around the arteries inside of the arm. Mounted inside of the cuff is a special gauge that is capable of sensing the small *oscillations* (fluctuations) in pressure. These fluctuations are produced as a result of the pressure the cuff has placed around the arm and are, in a basic sense, the arteries contracting (getting larger and smaller) with each heart beat. The monitor then measures how high and how quickly these contractions occur and converts that information into a digital value. A special valve, called an *electrode*, deflates the monitor automatically when the measurement is complete.

## Main Unit and Accessories

### Main Unit

### Accessories

4 AA Alkaline Batteries

## Special Features

### Noise Control System

A unique low-noise, rolling motor pump allows our monitors to offer quiet inflation.

### Large Three Row Display

A large LCD with three rows displaying time, date, and your blood pressure - systolic over diastolic - makes measuring your blood pressure easy to read.

### Memory Feature

Healthy Living® special memory feature allows the monitor to store up to 60 measurements along with the time and date for quick and easy comparisons.

### Time and Date Indications

This monitor automatically keeps track of the time, month and day.


## Accuracy

Healthy Living<sup>®</sup> monitors have been clinically tested against a scientific device called a *sphygmomanometer*, considered the gold standard in blood pressure measurement. All Healthy Living<sup>®</sup> monitors have performed equivalent to measurements taken with this scientific device and are within the accuracy limits prescribed by the American National Standard for Electronic or Automated Sphygmomanometers.

## Special Conditions

All Healthy Living<sup>®</sup> blood pressure monitors are safe and accurate to use with arrhythmias such as atrial or ventricular premature beats and atrial fibrillation.

## Digital Display Symbols


## Battery Installation

When the display shows , it's time to replace the batteries. This monitor uses 4 AA batteries.

Insert the batteries into the monitor making certain to align the polarities correctly.

Close the battery cover securely. It will snap into place when it is fully closed.


\* It's a good idea to remove the batteries if you know you will not be using your monitor for an extended amount of time.

## AC Capability

This Healthy Living® blood pressure monitor is equipped with AC capability should you choose to purchase an AC adapter to power your unit instead of batteries.

AC Adapter Output: 6V DC 1A

\*It is suggested that you remove the batteries from the compartment if you are going to use AC adapter for an extended amount of time.

## Selecting Time and Date

Your Healthy Living® monitor will keep track of the time, month, and day unless the batteries have run down or are removed from the monitor.


To set the time simply press the "MODE" button. The month will begin to flash. Press the "SET" button until the correct month appears on the display, then press the "MODE" button to store, follow this same process to set the day and time.

## Important Reminders

Follow your doctor's advice on how and when to take your blood pressure. Otherwise, here are some important hints:

### ***Before Taking a Measurement***

- ✓ Sit quietly for about 10 minutes. This will allow your body to return to its normal, resting state.
- ✓ Take your blood pressure in a comfortable environment - not too cold, not too hot - as temperature may affect your measurement.
- ✓ If you've just consumed beverages containing caffeine such as coffee, non-herbal tea, or cola, wait 30-45 minutes. Also don't smoke just before taking a measurement.


- ✓ Plug the cuff tube into the monitor before taking a measurement.

## Important Reminders



### *While Taking a Measurement*

- ✓ Sit still and quietly while measuring. Talking or moving may elevate measurements.
- ✓ For consistency, it's a good idea to measure your blood pressure using the same arm and at about the same time each day. If possible, we suggest using your left arm.
- ✓ When taking multiple measurements right after each other, make sure you wait at least 10 minutes in between. Waiting allows your blood vessels to return to their normal state.
- ✓ Sit with your legs uncrossed and your feet flat on the floor. Do not touch the cuff or monitor at any time during the measurement. Relax.

## Applying the Cuff

- Before you put the cuff on your arm make certain that you are not wearing a thick coat or sweater.
- If you are wearing a long sleeved shirt, roll up the sleeve, but make sure that it is not constricting circulation to your arm.
- The cuff should be positioned on your bicep and the red arterial marker should be in the center just above your elbow. Instructions are also printed directly on the cuff for reference.
- Wrap the cuff around your upper arm (usually the left) so it fits securely but is not too tight. You should be able to place 2 fingers in between the cuff and your arm. The velcro tab should fall within the "OK" range printed on the cuff for proper use.
- Place your arm on a table with your palm facing upwards. The inflation tube should line up with your arm. Do not allow the tube to kink or wrap around anything; it should fall down the center of your arm.

## Taking a Blood Pressure Measurement

- Once the cuff is properly applied and you are seated with your feet flat on the floor you are ready to take a measurement.

- Press the "START/STOP" button. The monitor will inflate to 160mmHg. If the monitor senses that you require a higher inflation level, it will automatically increase the air pressure to the level that's right for you.

- You will feel the cuff get tighter around your arm - this is normal. It should not inflate to the point of discomfort. If this occurs, the cuff is secured too tightly. Press "START/STOP" and rewrap the cuff leaving more room between your arm and the cuff.
- Numbers will appear on the digital display. This is the normal process the monitor goes through to determine your blood pressure. Don't be alarmed if you notice the numbers quickly change, or go higher than you believe your blood pressure to be.

## Taking a Blood Pressure Measurement

- When your measurement is complete the monitor will deflate automatically, releasing the air pressure in the cuff.

- The entire measurement process should take less than one minute. Once measurement is complete the display will flash between your blood pressure and pulse measurements.

**NOTE:** Blood pressure constantly changes during the course of a day. Occasionally, high or low readings do not necessarily mean your monitor is malfunctioning or that you are ill. It is important to speak with your doctor to learn what your normal blood pressure is and if you experience related abnormal measurements.

## Sixty Measurement Memory

When you take a measurement the monitor automatically stores it in memory space number one. As you continue to take measurements the monitor will push the previous measurement back one space. In other words, the value in space one will always be the most recent blood pressure measurement taken. Your monitor will store up to 60 measurements. To review measurements in the memory, press the "MEMORY" button. Press the "MEMORY" button in succession to review all 60 of your last blood pressure measurements.

## Taking Care of Your Blood Pressure Monitor

It's important to treat your monitor with care. Here are some basic tips:

- ✓ Keep the monitor out of direct sunlight, high humidity, and extreme temperature changes. 
- ✓ Try not to drop it on the floor.
- ✓ Do not take the monitor or the cuff apart. None of the monitor's parts are fixable at home. If your monitor needs repairs see details in warranty section.
- ✓ Keep the monitor free of dust and other materials. Avoid getting any part of the monitor or cuff extremely wet. If the monitor needs cleaning it's best to wipe it with a cloth slightly dampened with water. Healthy Living® recommends cleaning your monitor about once a month.
- ✓ Do not harshly fold or bend the cuff as this may damage it on the inside. 
- ✓ Do not inflate the monitor's cuff if it is not wrapped around your arm.
- ✓ Always remove the cuff tube from the monitor with care.
- ✓ For your safety, do not wrap the cuff around any other body part other than your arm.

## Troubleshooting

Error/Symbol	What Does It Mean?	How To Fix It
<b>ERROR</b>	Monitor may have been inflated without the cuff around the arm <i>OR</i> Cuff may not be positioned correctly on the arm.	Read page 17 (Applying the Cuff).
	You may have moved while the cuff was inflating, it is now attempting to reinflate. Check the cuff to make sure it's positioned correctly.	Read pages 15 & 16 (Important Reminders).
	You may have moved or spoken during the measurement process.	Remember that it is important not to move and speak as it may elevate measurements.
	You may have been startled by a noise or vibration while measuring.	Make sure you are in a quiet and comfortable place when taking a measurement.

## Troubleshooting

Error/Symbol	What Does It Mean?	How To Fix It
<b>ERROR</b>	You have been ignoring the flashing replace battery symbol (  ).	Replace with four new "AA" batteries. Make sure the polarities are positioned correctly in the battery compartment before continuing to measure your blood pressure.
	The batteries have run down.	Replace with 4 AA batteries quickly so you may continue to measure your blood pressure without any interruptions.
Blood pressure value seems too high or too low	Cuff may not be aligned correctly.	Read page 17 (Applying the Cuff).
The monitor keeps reinflating	Circuitry locked.	Remove then reinsert the batteries and then proceed to take measurement again.
Unable to take a measurement after receiving an error symbol.	The cuff may not be positioned correctly.	Read page 17 (Applying the Cuff).

## Five Year Limited Warranty

Samsung America, Inc. will, at its option, repair or replace your monitor, including the cuff, with the same or comparable model free of charge (except you must pay for shipping charges set forth below) for a period of (5) five years from the date of original purchase in the event of a defect in materials or workmanship occurring during normal use. This product is for home use only and is not intended for commercial or clinical use. This warranty only extends to the original retail purchaser. Purchase receipt or other proof of date of original purchase is required before warranty performance. The warranty on any replacement product provided under this warranty shall be for the unexpired portion of the warranty period applicable to the original product.

This warranty extends solely to failures due to defects in materials or workmanship occurring during normal use. It does not cover normal wear of the product. This warranty is void if: the product housing has been removed, if the product's label, logo or serial number have been removed or if the product fails to function properly as a result of accident, misuse, abuse, neglect, mishandling, misapplication, defective batteries, faulty installation, setup, adjustments, improper maintenance, alteration, maladjustment of controls, modification, power surges, commercial or clinical use of the product, service by anyone other than an authorized service center or acts beyond the manufacturer's or distributor's control.

The warranty and remedy provided herein are exclusive and in lieu of all other express and implied warranties and unless stated herein, any statements or representations made by any other person or entity are void. NEITHER THE MANUFACTURER, DISTRIBUTOR NOR THEIR AFFILIATES SHALL BE LIABLE FOR ANY LOSS, INCONVENIENCE, INJURY, OR DAMAGE INCLUDING DIRECT, SPECIAL, INCIDENTAL, OR CONSEQUENTIAL DAMAGES RESULTING FROM THE USE, IMPROPER USE, FAILURE TO HEED ANY WARNING, DISCLAIMER

## Five Year Limited Warranty

OR INSTRUCTION OR INABILITY TO USE THIS PRODUCT, WHETHER ARISING OUT OF ANY BREACH OF THIS WARRANTY OR ANY OTHER LEGAL THEORY. ALL EXPRESS AND IMPLIED WARRANTIES, INCLUDING THE WARRANTIES OF MERCHANTABILITY, ARE LIMITED TO THE APPLICABLE WARRANTY PERIOD SET FORTH ABOVE. Manufacturer and Distributor assume no responsibility for injuries, damages or penalties incurred or suffered resulting from the use of this product in an improper manner or manner or location other than for that for which it is intended or from failure to heed any warning, disclaimer or instruction.

Some states do not allow the exclusion or limitation of incidental or consequential damages, or limitations on how long an implied warranty lasts, so the above exclusions or limitations may not apply to you. This warranty gives your specific legal rights and you may also have other rights which vary from state to state.

Should your product prove defective during this warranty period, please return the product in the original carton prepaid to:

**SAMSUNG AMERICA, INC.**  
**c/o Healthy Living Division**  
**105 Challenger Road, 6th Floor**  
**Ridgefield Park, NJ 07660**  
**Attn: Customer Service**

Include \$11.95 for return shipping & insurance. (CHECK OR MONEY ORDER ONLY PLEASE, NON-REFUNDABLE). Sorry, we cannot deliver to P.O. Boxes. Please allow 6-8 weeks for delivery.

Include with your monitor a letter stating: your name, street address, daytime phone number, detailed description of the problem, a copy of your dated cash register receipt or other proof of purchase.

## Specifications

<b>Model</b>	<b>BA-507S</b>
Monitor Type	Automatic Inflate
Measurement Method	Oscillometric
Display	Liquid crystal digital display
Measurement Range	Pressure: 0-260 mm Hg Pulse: 40-199 beats/Minute
Accuracy	Pressure: within ± 3mm Hg Pulse: within ± 5%
Memory	Stores up to 60 measurements
Automatic Power Off	About 1 minute after measurement is complete
Battery Type	4 AA batteries (Alkaline)
Battery Life	About 6 months when used once a day, less when using a large cuff
Storage Temperature	Temperature: -4°F to 140°F Relative Humidity: less than 95%
Operating Temperature	Temperature: 50°F to 104°F Relative Humidity: less than 95%
Dimensions (approx.)	Monitor: 4 1/4" (w) x 5 3/4" (l) x 3" (h) Cuff: fits arms 9"-13" around
Weight (approx.)	14 oz (including batteries)

***For purposes of improvement, specifications are subject to change without notice.***

## Blood Pressure Diary

Date	Time	Systolic/Diastolic	Pulse
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	

## Blood Pressure Diary

**Date**      **Time**      **Systolic/Diastolic**      **Pulse**

		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	

## Blood Pressure Diary

**Date**      **Time**      **Systolic/Diastolic**      **Pulse**

		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	

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