Manual Inflate Blood Pressure Monitor

- Manually Controlled Inflation
- Handy, Lightweight Design
- 99 Memory Total
- Memory Average of All Stored Measurements
- Time and Date Indications
- Pulse Display

Important Warnings/Medical Disclaimer

Before using your blood pressure monitor, you should read and understand all instructions and follow all warnings.

Information in this manual is provided for informational purposes only. This product and manual are not meant to be a substitute for the advice provided by your own physician or other medical professional. You should not use this product, or the information contained herein, for diagnosing or treating a health problem, disease, or prescribing any medication. If you have or suspect that you have a medical problem, promptly contact your medical provider.

Table of Contents

Important Warnings/Medical Disclaimer _____	1
Healthy Living® _____	3-4
Commonly Asked Questions About Blood Pressure _____	5-8
Intended Use _____	8
Measurement Method _____	9
Important information before using the unit _____	10-11
Main Unit and Accessories _____	12
Special Features _____	13
Accuracy _____	13
Special Conditions _____	14
Battery Installation _____	14
Digital Display Symbols _____	15
Setting Date and Time _____	16
Applying the Cuff _____	17
Taking a Blood Pressure Measurement _____	18-19
Storing, recalling, and erasing measurement data _____	20-21
Memory - 99 Readings _____	21
Taking Care of Your Blood Pressure Monitor _____	22-23
Troubleshooting _____	24-25
5 Year Limited Warranty _____	26-27
Specifications _____	28
Blood Pressure Diary _____	29-32

Healthy Living®

Purchasing this blood pressure monitor is your first step in living a healthier life, but it is only the beginning. It is also important to:

STAY ACTIVE - Start out by taking just a short walk each day. Over time you'll begin to feel more energetic. Walking is a fun, easy way to help lose weight and also aids in building stronger bones and a healthier heart!

EAT HEALTHY - Adding more fruits and vegetables to your diet is a great way to start. Also, look at incorporating more whole grains, fish and low fat dairy products and keep a check on your sodium intake. If you're using a blood pressure monitor, you should already be aware of how sodium can affect your blood pressure. As always, consult your physician before beginning any dietary changes.

STAY POSITIVE - A healthy mind goes hand in hand with a healthy body. Daily relaxation techniques such as meditation and yoga may be just what you need to alleviate stress. Keep your mind focused on the positive things in life. Start out by giving yourself a pat on the back for those walks you are now taking or that apple you just ate!

WATCH YOUR WEIGHT - Sometimes just a small amount of weight loss may have a positive effect on your blood pressure. If you have already started exercising and eating healthier, the weight loss may just happen on its own!

REMEMBER: Consult your physician before you start or change your exercise routine, diet, vitamin/supplement intake and/or blood pressure monitoring.

Commonly Asked Questions About Blood Pressure

Why Should I Monitor My Blood Pressure At Home?

By monitoring at home you can manage your health care between visits to the doctor. Also, you now have the means to keep a record of your blood pressure measurements, a valuable tool for both you and your healthcare provider.


Home monitoring may help alleviate White Coat Hypertension. This term refers to what may cause a person's blood pressure to be higher when taken in a medical setting. At times the stress of just being in a doctor's office may cause blood pressure to rise. When monitoring at home you may get a truer indication of your blood pressure because you are in a relaxed, comfortable space.

Warning: Do not use your results from your home monitoring for self-diagnosis or treatment. Please consult your medical provider about your health or condition.

What is Blood Pressure?


The term "blood pressure" refers to the force or pressure that is exerted on the body's arteries as blood flows through them. Each time the heart beats or contracts, it produces pressure in the arteries. This blood pressure is what moves blood through the body, supplying oxygen and nutrients to every organ.

Commonly Asked Questions About Blood Pressure

When the heart contracts it sends blood out into the body. This pressure is referred to as **systolic**.

The pressure that is created as the heart pulls blood back through is called **diastolic**.


Blood Pressure is recorded as **systolic over diastolic**.

120 80

What is considered Normal Blood Pressure?

The Seventh Report of the Joint National Committee on Prevention, Detection, Evaluation, and Treatment of High Blood Pressure from the National Heart, Lung and Blood Institute has set forth the following guidelines regarding high blood pressure. Guidelines provided are for information purposes only. Please consult with your physician for proper diagnosis.

Commonly Asked Questions About Blood Pressure

Why Does My Blood Pressure Fluctuate?

Blood pressure changes with every beat of the heart and is in constant fluctuation throughout the day. In addition to natural changes, other factors that may affect blood pressure are:

- Temperature (too hot or cold)
- Humidity
- Restlessness
- Lack of sleep
- Physical exercise
- Constipation
- Fear
- Anger
- Anxiety
- Food / beverage consumption
- Your physical condition
- Stress

Be sure to speak with your healthcare professional to determine what is affecting your blood pressure.

How Do I Know If I Have High Blood Pressure?

People have the misconception that they can "feel" their blood pressure rise. In reality it is nearly impossible to tell if your pressure is elevated without having it checked. This is why hypertension has been named "the silent killer". Start monitoring your blood pressure today and discuss the results with your healthcare professional.

Commonly Asked Questions About Blood Pressure

Why Is It Important To Control My Blood Pressure?

Left uncontrolled hypertension can place an incredible amount of stress on your heart and arteries. They are both forced to work harder to keep blood flowing through your body. Over time this stress may result in health related problems including enlargement of the heart, loss of artery elasticity, heart disease, stroke or death. Making simple lifestyle changes now may help you stay healthy in the future.

Intended Use


Healthy Living® blood pressure monitors are designed for home use. They are noninvasive; no part enters your body, and provide systolic, diastolic and pulse rate measurements. These monitors are not recommended for use by people under 18 years of age and are not suitable for clinical use. This monitor is not intended to be a diagnostic device.

Measurement Method

Healthy Living[®] monitors measure blood pressure and pulse rate by using an oscillometric method, meaning the fluctuations in pressure are measured. Once the cuff is wrapped around the upper arm, the monitor is turned on and you begin to inflate the cuff. The inflation of the cuff creates pressure around the arteries inside the arm. Within the cuff is a gauge which senses the fluctuations (oscillations) in pressure. The fluctuations measured are the arteries contracting with each heart beat as a result of the pressure the cuff has placed on the arm. The monitor measures these contractions and converts the information to a digital value. That is the result displayed on the monitor screen.


Important information before using the unit

- Blood pressure measurements should be interpreted by a physician or trained health professional who is familiar with one's medical history. By using the unit regularly and recording the results for the physician to interpret, one can keep the physician informed of the continuing trends in one's blood pressure.
- Sit quietly for about 10 minutes. This will allow your body to return to its normal resting state.
- Wrap the cuff snug around your upper arm and lay your arm on a flat surface at the same level as your heart. Your feet should be flat on the floor.
- The pumping of the bulb should be the only movement you make during the measurement. Pump with consistent, steady pressure.
- Perform the measurement in a quiet environment, in a relaxed position at room temperature. For consistency, we recommend using the same arm (left arm if possible) and testing around the same time each day.

Important information before using the unit


- If taking multiple measurements, please allow 10 minutes between measurements. This will allow your blood vessels to return to their normal state.
- If you have consumed beverages containing caffeine or have smoked, please allow 30 - 45 minutes before taking your measurement.
- Do not place the cuff over a jacket or sweater sleeve; it must be wrapped around your exposed upper arm.
- Keep in mind that blood pressure naturally varies from time to time throughout the day and is affected by many different factors, such as smoking, alcohol consumption, medication, and physical activity.
- Make sure the cuff tube is attached to the connector just above the bulb as shown above.

Main Unit and Accessories

Main Unit

(See page 15 for descriptions of screen icons.)

Accessories

- One 3V CR2032 Battery
- 9" ~ 13" Regular Cuff

Special Features

Heart Sense®

The Heart Sense® Inflation System senses your blood pressure before the measurement is complete, signaling the monitor to adjust the air pressure amount. The result is a quick, comfortable measurement each time you check your blood pressure.

Easy to Hold Monitor

This compact monitor puts tracking your health in your hands.

Memory Feature

Your Healthy Living® Manual Inflate Monitor holds up to 99 readings and will also average all the readings in memory!

Time and Date

The time and date will appear as long as the batteries are in or an AC Adapter is used. Once the batteries are removed or the adapter unplugged the time and date will have to be reset. The data in memory will remain.

Accuracy

Healthy Living® monitors have been clinically tested against a scientific device called a *sphygmomanometer*, considered the gold standard in blood pressure measurement. All Healthy Living® monitors have performed equivalent to measurements taken with this scientific device and are within the accuracy limits prescribed by the American National Standard for Electronic or Automated Sphygmomanometers. Healthy Living® monitors have been clinically tested and have met all the requirements to achieve AAMI (Association for the Advancement of Medical Instrumentation) certification.



Special Conditions

All Healthy Living® BPMs are safe and accurate to use with arrhythmias such as atrial or ventricular premature beats and atrial fibrillation. We recommend, however, that if you have these conditions or any other special conditions you should check with your physician before using your Healthy Living® BPM.

Battery Installation

Locate the battery cover on the backside of the monitor face. Using a coin or similar object, twist the latch in the direction of the arrow. Once the latch is off the battery should easily flip out. To replace, simply drop in a new 3V CR2032 battery making sure the positive side is facing up. The positive side will have a "+" stamped on it. Place the latch back in place and twist in the opposite direction of the arrow to tighten.

When the display shows "  " it is time to replace with a new battery.


* To make the batteries last longer, you may want to remove them if you are not using your monitor for an extended period of time.

* Keep batteries out of reach of children.

Digital Display Symbols

Measurement Reading Screen


Setting Date and Time

- a. Press "MODE" button ("month" starts flashing)
Press the set button to set and adjust the correct month (1, 2, 3 ...12)
- b. Press "MODE" button again ("date" starts flashing)
Press the set button to set and adjust the correct date
- c. Press "MODE" button again ("hour" starts flashing)
Press the set button button to set and adjust the correct hour in 12-hour format. A "PM" will appear for PM and disappear for AM.
- d. Press "MODE" button again ("minute" starts flashing)
Press the set button button to set and adjust the correct minute (00, 01, 02 ...59)

After entering these settings you will be returned to the Rest mode and can take a measurement by pressing START/STOP.

Applying the Cuff

- Insert the Cuff Plug into the Cuff Plug Jack on the Monitor.
- Pull the Velcro tab so the cuff is laying flat with the Velcro facing up at you.
- Pull the cuff to the left of the metal ring to form an opening for your arm.
- Put your left arm through the cuff loop. The bottom of the cuff should be approx. 0.5 inch above your elbow. The cuff tube should run along the center of your arm towards your upturned palm.
- If the cuff is properly positioned the tube should be running down your arm as in the illustration to the right.
- Pull the cuff so that the top and bottom edges are snug around your arm. Do not tighten to the point of being uncomfortable.
- When the cuff is positioned properly, press the Velcro firmly against the pile side of the cuff.
- Sit in a chair and place your arm on the table so that the cuff is at the same level as your heart.
- Relax your arm and sit with your palm facing upward.
- Assure there are no kinks in the air tube.


Taking a Blood Pressure Measurement

1. Once the cuff is properly applied and you are seated with your feet flat on the floor, you are ready to take a measurement.

2. Holding the pump bulb in your hand, press the START/STOP button, the "🌀" icon will appear. This tells you to release any trapped air in the cuff by pressing and holding the air release knob just above the bulb.

* **ERROR** will appear and 4 beeps will sound if the air is not released before pumping. The cuff must be fully deflated.

3. When the "🌀" icon begins flashing start squeezing the Pump Bulb. This will start to inflate the cuff and you will feel the cuff get tighter.

4. The cuff should not inflate to the point of discomfort, if this happens the cuff is wrapped too tightly. Depress the air release knob at the base of the Pump Bulb to release the air and press the "START/STOP" button to stop the measurement. Rewrap the cuff leaving more room between your arm and the cuff.

5. Continue to squeeze the Pump Bulb until you hear a beep. This is your signal to STOP squeezing and the "🌀" icon will be displayed, sit still and quiet while measuring. A flashing "♥" appears as your measurement is being taken and the bulb symbol will stop flashing. **(Do not inflate the monitor higher than 240 mmHg).**

Taking a Blood Pressure Measurement

*If the bulb symbol again begins to flash, that is an indication that more air is needed and you should resume pumping until the flashing stops.

6. As the monitor determines your blood pressure, numbers will appear and change on the display screen until your reading is complete.


Inflation

7. When your measurement is complete the " " icon will appear. This is your signal to deflate the cuff by depressing the Air Release knob.

Deflation

8. The process should take less than a minute and once the measurement is complete your blood pressure and pulse rate will be displayed. The monitor will turn itself off approximately 1 minute after completing a reading.

Storing, recalling, and erasing measurement data

• Storing data

*After each measurement, the systolic and diastolic rates, pulse rate and the time & date will be automatically stored. The memory holds data for the latest 99 measurements. If more than 99 measurements are entered, the memory will automatically clear out the oldest data.

• Recalling data

*Press the "MEMORY" button to get your stored memory data. If there is no data in the memory nothing will appear.

*If there is data in memory, the first reading will be an average of all the readings in memory (systolic and diastolic pressure and pulse rate).

*Press the "MEMORY" button again and the latest recorded set of measurements will appear along with the recording time and date.

*Press the "Memory" button to read the next set of stored data. The data displayed starts with the most recent reading and continues backwards to the oldest.

Storing, recalling, and erasing measurement data

• Erasing data

*Press "MEMORY" into Memory mode.

*Press and hold on "MODE" button and then press "SET" button, all measurement data will be erased.

*To confirm the data has been erased, press the "MEMORY" button and no data should appear.

Memory - 99 Readings

When you take a measurement, the monitor automatically stores it in memory space number one. As you continue to take measurements, the monitor will push the previous measurement back one space. In other words, the value in space one will always be the most recent blood pressure measurement taken.

Taking Care of Your Blood Pressure Monitor

Maintenance

Battery Removal

If the battery has been removed from your monitor, your stored data will not be affected. You will have to reset the time and date once power is restored, but the data in the memory zones will remain.

Dispose of batteries in accordance with applicable local law.

Use a damp cloth to wipe the monitor and dry immediately with a dry cloth.

Use only a dry cloth to wipe the cuff.

Do not use strong cleaning agents for cleaning.

When the unit is not to be used for a long time, it is best to remove the battery.

Safe Keeping

Do not store the monitor in direct sunlight, high temperatures, or in humid or dusty places.

Do not store in extremely low (less than -4 °F / -20 °C) or high (more than 122 °F / 50 °C) temperatures.

Taking Care of Your Blood Pressure Monitor

To prevent malfunctions

Use the cuff only as instructed and be careful not to twist it. Do not inflate the cuff when it is not wrapped around the arm.

Do not attempt to disassemble or change any monitor parts, including the cuff.

Treat the monitor with care, do not drop it and prevent all strong impacts.

Warnings

The cuff size is suitable for adult arms 9 ~ 13 inches around. The device is not to be used if your arm has a wound or injury.

In the case that the cuff does not stop inflating, remove the cuff at once.

Do not disassemble or modify the device in any way. This could result in inaccurate readings and voids the device warranty.


Do not use the cuff on any body part except your upper arm. Keep out of reach of children.

Do not use your monitor to diagnose or treat any health problems or illnesses, or prescribe any medication. Consult your doctor about your health or condition.

Troubleshooting

Error/Symbol	What Does It Mean?	How To Fix It
<p>ERROR</p>	<p>Monitor may have been inflated without the cuff around the arm <i>OR</i> Cuff may not be positioned correctly on the arm.</p>	<p>Read page 17 (Applying the Cuff).</p>
	<p>The cuff tube may be folded <i>OR</i> You may have moved while the cuff was inflating.</p>	<p>Read pages 10 & 11 (Important information before using the unit).</p>
	<p>You may have moved or spoken during the measurement process.</p>	<p>Remember that it is important not to move and speak as it may elevate measurements.</p>
	<p>You may have been startled by a noise or vibration while measuring.</p>	<p>Make sure you are in a quiet and comfortable place when taking a measurement.</p>
<p>ERROR</p>	<p>All air has not been released from the cuff, it is not fully deflated.</p>	<p>Deflate the monitor by depressing the Air Release Knob.</p>
	<p>The monitor is inflated higher than 300 mmHg</p>	
<p>continued beep sound</p>		

Troubleshooting

Error/Symbol	What Does It Mean?	How To Fix It
	<p>The batteries have run down.</p>	<p>Replace with 1 CR2032 battery.</p>
<p>Blood pressure value seems too high or too low</p>	<p>Cuff may not be aligned correctly</p>	<p>Read page 17 (Applying the Cuff).</p>
	<p>Circuitry locked.</p>	<p>Remove and reinsert the batteries and then proceed to take measurement again.</p>
<p>Unable to take a measurement after receiving an error symbol.</p>	<p>The cuff may not be positioned correctly.</p>	<p>Read page 17 (Applying the Cuff).</p>
	<p>EP (Electronic Prom) error may caused by calibration setting and memory problem.</p>	<p>Please call Customer Service for assistance.</p>

5 Year Limited Warranty

This Blood Pressure Monitor is warranted to be free from defects in material and workmanship for five (5) years from the date of the original purchase under normal household use. In the event of a defect in materials or workmanship occurring during normal household use, Samsung America, Inc. ("Samsung America") will, at its option, repair or replace your monitor, including the cuff, with the same or comparable model free of charge (except you must pay for shipping charges) for five (5) years from the date of purchase. This product is for home use only and is not intended for commercial or clinical use. This warranty only extends to the original retail purchaser. Purchase receipt or other proof of original purchase is required before warranty performance. The warranty period on any replacement or repaired product provided under this warranty shall be for the unexpired portion of the warranty period applicable to the original product.

This warranty extends solely to failures due to defects in materials or workmanship occurring during normal use. It does not cover normal wear of the product. This warranty is void if: the product housing has been removed, if the product's label, logo or serial number have been removed or if the product fails to function properly as a result of accident, misuse, abuse, neglect, mishandling, misapplication, defective batteries, faulty installation, setup, adjustments, improper maintenance, alteration, maladjustment of controls, modification, power surges, commercial or clinical use of the product, improper or incorrectly performed service, or acts beyond the manufacturer's or Samsung America's control.

The warranty and remedy provided herein are exclusive and in lieu of all other express and implied warranties and unless stated herein, any statements or representations made by any other person or entity are void. NEITHER THE MANUFACTURER, SAMSUNG AMERICA, NOR THEIR AFFILIATES SHALL BE LIABLE FOR ANY LOSS, INCONVENIENCE, INJURY, OR DAMAGE INCLUDING DIRECT, SPECIAL, INCIDENTAL, OR CONSEQUENTIAL DAMAGES RESULTING FROM THE USE, IMPROPER USE, FAILURE TO HEED ANY WARNING, DISCLAIMER OR INSTRUCTION OR INABILITY TO USE THIS PRODUCT, WHETHER

5 Year Limited Warranty

ARISING OUT OF ANY BREACH OF THIS WARRANTY OR ANY OTHER LEGAL THEORY. ALL EXPRESS AND IMPLIED WARRANTIES, INCLUDING THE WARRANTIES OF MERCHANTABILITY, ARE LIMITED TO THE APPLICABLE WARRANTY PERIOD SET FORTH ABOVE. Manufacturer and Samsung America assume no responsibility for injuries, damages or penalties incurred or suffered resulting from the use of this product in an improper manner or location other than for that for which it is intended, or from failure to heed any warning, disclaimer or instruction.

Some states do not allow the exclusion or limitation of incidental or consequential damages, or limitations on how long an implied warranty lasts, so the above exclusions or limitations may not apply to you. This warranty gives you specific legal rights and you may also have other rights which vary from state to state.

Should your product prove defective during this warranty period, please return the product in the original carton prepaid to:

SAMSUNG AMERICA, INC.
c/o Healthy Living Division
105 Challenger Road, 6th Floor
Ridgefield Park, NJ 07660
Attn: Customer Service

Include \$11.95 for return shipping & insurance within the continental U.S. (CHECK OR MONEY ORDER ONLY PLEASE, NON-REFUNDABLE). **If you are not the original purchaser of the unit with proof of purchase, or if your warranty period has expired, the price for repairing or replacing may vary.** Sorry, we cannot deliver to P.O. Boxes *or outside the U.S.*

Please allow 6-8 weeks for delivery.

Include with your monitor a letter stating: your name, street address, day-time phone number, detailed description of the problem, a copy of your dated cash register receipt or other proof of purchase.

This warranty applies only to products purchased in the U.S. For products purchased in, but used outside, the U.S., this warranty covers only warranty services within the U.S. (and does not include shipping outside the U.S.).

Specifications

Model	BHM-1008
Monitor	Manual Inflate
Measurement Method	Oscillometric
Display	Liquid crystal digital display
Measurement Range	Pressure: 0 ~ 300 mmHg Pulse: 40 ~ 199 beats/Minute
Accuracy	Pressure: within ± 3 mmHg Pulse: within $\pm 5\%$
Memory	Stores up to 99 measurements
Automatic Power Off	Approximately 1 minute after measurement is complete
Battery Type	One 3V CR2032
Battery Life	Approximately 180 times measurement (once a day measurement)
Storage Temperature	Temperature: -4 °F ~ 122 °F/-20 °C ~ 50 °C Relative Humidity: 10 ~ 85% R.H.
Operating Temperature	Temperature: 50 °F ~ 104 °F/10 °C ~ 40 °C Relative Humidity: 30 ~ 85% R.H.
Dimensions (approx.)	Monitor: 6 $\frac{13}{16}$ " (L) x 2 $\frac{13}{16}$ " (W) x 1 $\frac{3}{4}$ " (H) Standard Cuff: fits arms 9" - 13" around
Weight (approx.)	3.53 oz (including battery)

For purposes of improvement, specifications are subject to change without notice.

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Copyright © 2006 Samsung America, Inc. Todos los derechos reservados.
Copyright © 2007 Samsung America, Inc.